# STORZ

## KARL STORZ—ENDOSKOPE

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**Important information  
for users of KARL STORZ  
instruments**

It is recommended to check the suitability of the product for the intended procedure prior to use.

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.


This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the ThermoFlator®. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment.  
Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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②: The KARL STORZ-SCB interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

① **Note:** If two or more identical devices are used in the KARL STORZ-SCB system, please observe the section "Identical Devices on SCB" on page 30.

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**Controls, displays,  
connectors, and  
their uses**

- ① Main power switch
- ② Gas supply bar graph
- ③ Patient line pressure actual-value display
- ④ Patient line pressure setpoint display
- ⑤ Digital setpoint/actual-value display of patient line pressure
- ⑥ Patient pressure mode indicator
- ⑦ Gas flow actual-value display
- ⑧ Gas flow setpoint display
- ⑨ Digital gas flow setpoint/actual-value display
- ⑩ Gas volume display
- ⑪ Blow-off valve indicator
- ⑫ Heating element indicator
- ⑬ Insufflation line fitting, line to patient
- ⑭ Optitherm® heating element connector
- ⑮ Volume reset button
- ⑯ Insufflation START/STOP button
- ⑰ Flow mode indicator
- ⑱ Flow setting switches
- ⑲ Flow mode switch (M button)
- ㉐ Drawer for operator's quick reference guide
- ㉑ Pressure setting switches
- ㉒ Insufflation initialization mode switch with indicator
- ㉓ Label strip (replaceable)
- ㉔ Optitherm® heating element
- ㉕ Gas inlet (American standard fitting)
- ㉖ CO<sub>2</sub> cylinder clamping fixture
- ㉗ SCB connectors
- ㉘ Potential equalization connector
- ㉙ Power fuses
- ㉚ Power cord socket

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## Symbols employed

	Follow instructions for use
	Gas cylinder pressure
	Patient pressure
	Insufflation pressure
	Gas Flow
	Semicontinuous flow mode
	Intermittent flow mode
	CO <sub>2</sub> Volume
	Blow-off valve indicator
	Heating element indicator
	Gas outlet to the patient
	Heating element connector
	Gas volume reset
	START/STOP
	Insufflation initialization mode
	Potential equalization connector
	Applied part of type BF
	Alternating current
	Electronic information product pollution control (China RoHS)
	This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Manufacturer

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## Contents

Images of the equipment .....	IV, VII
Controls, displays, connectors, and their uses .....	V
Symbols employed .....	VI
<b>General information</b>	
Device description .....	2
<b>Safety instructions</b>	
Warnings and cautions .....	3
Intended use .....	6
User qualification .....	8
Safety precautions at the site of installation .....	9
Safety precautions when operating the unit .....	9
Safety features .....	10
<b>Installation and operating instructions</b>	
Unpacking the equipment .....	13
Basic equipment .....	13
Installing and connecting up the unit .....	13
Connecting the CO <sub>2</sub> cylinder .....	15
Gas supply modes .....	16
Operating the unit .....	17
Heating up the CO <sub>2</sub> gas (Optitherm®) .....	20
Selecting the insufflation mode .....	22
Preselecting the gas flow .....	24
SECUVENT® blow-off valve .....	26
Test for proper operation .....	27
Initiating CO <sub>2</sub> insufflation .....	29
Service programs P1, P2 .....	30
<b>Maintenance</b>	
Fuse replacement .....	31
Cleaning, disinfection and sterilization .....	32
Maintenance and safety check .....	39
Servicing and repair .....	40
Repair program .....	41
Limitation of liability .....	42
Manufacturer's warranty .....	42
<b>Technical description</b>	
Troubleshooting .....	43
Technical data .....	44
Technical documentation .....	45
Insufflation operational schematic .....	46
<b>Spare parts, recommended accessories</b>	
Spare parts/accessories .....	47
Accessories .....	48
<b>Appendix</b>	
Electromagnetic Compatibility (EMC) Information .....	49

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### Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

### Device description

The ThermoFlator® is an insufflation device for universal application in laparoscopic, colonoscopic and thoracoscopic examinations and operations. State-of-the-art technology for the measurement and control of both pressure and flow enable the use of different operating modes, which can be tailored to specific situations such as the creation of pneumoperitoneum, the use of lasers or the performance of HF surgery. Easy handling and clearly arranged controls together with several safety circuits ensure maximum patient safety.

To compensate for the considerable loss of gas occurring in complex laparoscopic operations, the ThermoFlator® has been designed for a high flow rate of up to 30 l/min.

Use of the OPTITHERM® heating element can reduce the risk of hypothermia.

### Introduction of two different gas supply modes

The ThermoFlator® is frequently used on a central gas supply with pressures below 7 bar.

In older units, this used to trigger a warning because the pressure measuring system of the unit interpreted it as an ‘empty’ high-pressure bottle.

The new ThermoFlator® can distinguish between two different supply modes:

- high pressure and
- low pressure.

To set the desired service mode see section ‘Service programs’ on page 30.

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## Warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make these words stand out more clearly, they are accompanied by a pictogram.

**Warning:** A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

**Caution:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

**Note:** A Note indicates special information about operating the product, or clarifies important information.

**Warning:** Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

**Note:** Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

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## Safety instructions

### Warnings and cautions

**Warning:** Test this unit prior to each surgical procedure to ensure that it functions correctly.

**Warning:** DO NOT use the ThermoFlator® in the presence of flammable anesthetics.

**Warning:** To avoid the risk of electrical shock, do not open the unit. Refer servicing to qualified personnel.

**Warning:** Grounding reliability can only be achieved when the equipment is connected to 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

**Warning:** The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards.

**Warning:** Keep the unit out of reach of patients.

**Warning:** Always unplug the unit before all maintenance work (e.g. cleaning).

**Warning:** The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

**Warning:** Always keep a spare device to hand in case the first device should fail.

**Warning:** Make sure there are no external forces acting on the cavity. This can increase the intracavitary pressure and/or cause pressure fluctuations.

**Warning:** Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cf. IEC 60601-1-1).

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**Safety instructions**  
**Warnings and cautions**

- !** **Warning:** The device may only be operated under the ambient conditions specified.
- !** **Warning:** For safety reasons, do not simultaneously touch the device output sockets and the patient.
- !** **Caution:** Avoid allowing liquids to enter the unit. Do not store liquids on or above the unit.
- !** **Caution:** Only operate the unit at the line voltage range stated on the manufacturer's identification plate.
- !** **Caution:** When replacing fuses, use only fuses of the correct rating.
- !** **Caution:** To avoid contamination of the unit as a result of the reverse flow of gas or body fluid, a sterile CO<sub>2</sub> filter must be inserted between the insufflation tube and unit inlet.
- !** **Caution:** The CO<sub>2</sub> bottle must be placed vertically while operating the instrument, otherwise the functioning of the instrument cannot be guaranteed.
- !** **Caution:** CO<sub>2</sub> cylinders connected up to the unit must be secured to prevent them from falling over.
- i** **Note:** This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap. Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.
- i** **Note:** When disposing of accessories, the respective national regulations must be observed.

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## Intended use

CO<sub>2</sub> insufflators and their accessories are used to create and maintain the cavity during diagnostic endoscopic interventions as well as to expulse ambient air during open, endoscopically assisted or endoscopic interventions.

## Indication

The insufflator is used to create a cavity during the following diagnostic and therapeutic interventions:

- laparoscopy
- thoracoscopy
- endoscopy of the upper and lower gastrointestinal tract (e.g., TEO®, colonoscopy)

Only sterile CO<sub>2</sub> gas may be used.

Use of the device in fields other than those indicated above is not allowed for safety reasons.

The ThermoFlator® may only be operated with tubing sets and accessories which have been designated as suitable for the device by KARL STORZ.

**Caution:** If tube systems other than those prescribed are used, no guarantee can be given for safe operation of the instrument.

Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

**Warning:** additional devices which are connected to electronic medical devices must conform to the respective IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore, all configurations must satisfy the requirements for medical electrical systems (see IEC 60601-1-1 or paragraph 16 of the 3rd Ed. of IEC 60601-1). Any person who connects additional devices to electrical medical devices configures a medical system, and is therefore responsible for ensuring that this system satisfies the requirements for medical electrical systems. It must be taken into account that local regulations take priority over the aforementioned requirements. If in doubt, please contact your local representative or technical support.

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### **Contraindications for Use**

Contraindications arising directly from the patient findings must be observed.  
The unit is contraindicated for hysteroscopic applications.

### **CO<sub>2</sub> absorption**

CO<sub>2</sub> absorption occurs during insufflation. The body absorbs a portion of the CO<sub>2</sub> used for the insufflation. A too high CO<sub>2</sub> concentration in the blood or in the respiratory system can lead to the death of the patient in isolated cases. Observe vitals throughout the entire procedure and ensure sufficient patient respiration. High pressure or high gas flow encourage CO<sub>2</sub> absorption.

### **Changes to circulatory, metabolic and respiratory function**

The following symptoms can occur during CO<sub>2</sub> insufflation:

- Hypercapnia
- Acidosis
- Decrease in cardiac output
- Decrease in venous return

### **Embolism**

Incorrect positioning of the insufflation instrument or a high intraabdominal pressure can lead to a CO<sub>2</sub> embolism. Avoid high pressures and verify the correct positioning of the insufflation instrument.

### **Hypothermia**

The CO<sub>2</sub> insufflation can lead to a decrease in body temperature. Monitor the body temperature throughout the intervention.

### **Idiosyncratic reactions**

The risk of a metabolic imbalance due to increased CO<sub>2</sub> absorption is increased for patients with sickle cell anemia or pulmonary insufficiency.

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**Dehydration**

For longer operations with high gas levels, insufflation can lead to a drying out of the tissues or cause tissue damage. Unnecessary leaks must be avoided for this reason.

The decisive factor for the incidence of these complications is the length of the pneumoperitoneum, the level of the intraabdominal pressure and the CO<sub>2</sub> absorption associated with it. Increased pressure and flow values must thus be avoided. Pressures in excess of 15 mmHg are only required in isolated cases.

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**Precautions**

To prevent patient cross contamination, a hydrophobic bacterial filter must be used.

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
**User qualification**

The Thermoflator® may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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### **Safety precautions at the site of installation**

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies for easily combustible and explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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### **Safety precautions when operating the unit**

It is the user's responsibility to make sure the unit is safe and operates properly before using it (see the section 'Test for proper operation').

During treatment using the ThermoFlator®, the patient must be treated and kept under observation with the usual medical care. This includes keeping an endoscopic check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

Any treatment may only be performed if there is visual observation of the action of the unit.

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## Safety features

The ThermoFlator® has the following safety features:

- **Self-test procedure:** The self-test procedure is carried out every time the instrument is switched on. The following components are tested:
  - **General functional test:** A sequence of 6 short sounds are given to confirm correct operation.
  - **Checking pressure measurement:** If the result is negative a continuous warning sound is given and the message ‚Pr Err‘ appears in the displays ⑤/⑩.
  - **Checking flow measurement:** If the result is negative a continuous warning sound is given and the message ‚Fl Err‘ appears in the displays ⑨/⑩.
- **Note:** Alarms which sound when switching on the unit (self-test), can only be canceled by switching off the unit.
- **Checking the intracorporeal pressure**
  - **Overpressure:** As soon as the intracorporeal pressure exceeds the current setpoint by more than 5 mmHg for longer than 3 seconds, a continuous acoustic warning signal is given (alternating long and short sound) and the gas supply is stopped. After a further 1.5 seconds the overpressure is actively released by a blow-off valve. If the pressure exceeds the limit value of 33 mmHg for longer than 5 seconds a continuous alarm sound is given.
  - **Underpressure:** If an intracorporeal underpressure is detected a continuous acoustic warning signal (2 successive short warning sounds) is given.
  - **Dangerous pressure rise:** Three warning sounds in one-second rhythm are emitted if the pressure increase speed exceeds a certain value (2.5 mmHg/0.1 s) in the initialization mode. In this way attention is called to incorrect placement of the VERESS needle, e.g. in the tissue. Simultaneous to this acoustic warning signal, the digital pressure display ⑤ flashes.

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## Safety instructions

### Safety features

- **Gas flow interrupted:** If a flow = 0 is recognized during the operating mode with semi-continuous gas flow, a warning sound is given in one-second rhythm and the digital gas flow display ⑨ flashes. Since normally one always works with a gas leak in this operating mode, this warning can call the user's attention to a kinked insufflation tube or the like, and therefore to a false pressure indication.

- **Checking the gas supply** (high pressure mode) If the connected CO<sub>2</sub> cylinder is empty, a continuous acoustic warning signal is given as soon as insufflation is started (3 warning sounds in quick succession).

**i Note:** This does not apply to low-pressure gas supply mode.

- **Checking gas preheating** The internal gas preheating is monitored by a second circuit, independently of the control. If a fault arises the unit is switched off.

- **Checking the gas temperature to the patient**
  - **If the gas temperature is too high** (higher than 41 °C) at the outlet nozzle of the Optitherm® heating element a continuous alarm sounds and the temperature symbol ⑫ on the front panel lights up red. The heating element is switched off in this case. Insufflation with CO<sub>2</sub> is, however, still possible. The temperature alarm resets itself if the gas temperature falls again to values below 41 °C. Additionally, it can be shut off by pulling out the plug for the Optitherm® heating element.

- **If the temperature sensor of the** Optitherm® heating element is faulty a continuous alarm also sounds and the temperature symbol on the front panel lights up red.  
The device must be switched off and switched on again, either with a properly functioning Optitherm® heating element or with no heating element at all.

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**Safety instructions****Safety features**

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**• Passive safety measures**

- The pressure measurement system has a redundant stand-by.
- Both pressure control stages are equipped with safety pressure release valves.

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### Unpacking the equipment

Carefully unpack the ThermoFlator® and remove the unit and its accessories from their packaging.  
Check for missing items and evidence of shipping damage.  
File any complaints with the manufacturer or supplier immediately.  
If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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### Basic equipment

KARL STORZ Set No.: 264320 08-1

- 1 SCB® ThermoFlator® 264320 20-1
- 1 Power cord 400 A
- 1 Optitherm® heating element 204320 30
- 1 Silicone tubing set 204000 43, sterilizable
- 1 Universal wrench 204000 30
- 1 SCB connecting cord 200901 70, length 100 cm
- 1 Package of CO<sub>2</sub> filters, sterile, disposable, package of 10 pcs.
- 1 Instruction manual
- 1 Operator's quick reference guide (self-adhesive)

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
### Installing and connecting up the unit

**Note:** The ThermoFlator® including the accessories connected may be used only in medical facilities having electrical installations conforming to the applicable national, state, and local electrical codes.

**Warning:** The unit is not intended for use in potentially explosive atmospheres. Do not operate the unit within the demarcated hazard zone while explosive anesthetic gases are in use.

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### Fitting a label strip to the front panel

Undo the left-hand front side cover with a screwdriver.

**i Note:** The appropriate controls/indicators on the unit are designated on the label strip.

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### Potential equalization

**Warning:** Connect the potential equalization cable up to the potential equalization connector ⑧. The unit's ground line should be installed by a qualified electrician.

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Only operate the unit at the line voltage range stated on the manufacturer's identification plate.

Insert the power cord into the power socket ⑩ as far as it will go.

**Warning:** Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

**Note:** Please ensure that access to the power plug is ensured at the site of installation.

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### KARL STORZ-SCB

**Note:** To avoid mistakenly pulling out the SCB connecting cord, the SCB connector possesses a protection device.

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets ⑦. Connect the other end of the cord to the KARL STORZ-SCB (KARL STORZ Communication Bus) control unit or other SCB units (see KARL STORZ-SCB control NEO System Instruction Manual).

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### Attaching the CO<sub>2</sub> cylinder at the unit

Open the clamping fixture ⑥ located on the unit's rear panel by dropping its baseplate downward and then raise the upper clamping ring until it locks into place.

Insert the CO<sub>2</sub> cylinder in the clamping fixture vertically and close the locking clip (see safety instructions page 5).

**Warning:** Operate only with sterile CO<sub>2</sub> gas.

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### Connecting the CO<sub>2</sub> cylinder

Connect the CO<sub>2</sub> cylinder's pressure fitting to the unit's gas inlet ⑤ using the high-pressure tube provided.

**Note:** The unit's gas inlet ⑤ is a standard American type. Appropriate adaptor tubes will be needed to connect CO<sub>2</sub> cylinders having German-standard or keyed PIN-Index pressure fittings.

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### Connecting CO<sub>2</sub> cylinders equipped with German-standard/ISO fittings

Align the adapter tube's threaded ring on the CO<sub>2</sub> cylinder's pressure fittings and tighten down the ring using the universal wrench provided.

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### Connecting CO<sub>2</sub> cylinders equipped with keyed PIN-Index fittings

Align the adapter tube's threaded ring on the CO<sub>2</sub> cylinder's pressure fittings and tighten down the ring.

Use the wrench provided (see the illustration) to open the CO<sub>2</sub> cylinder's pressure valve.

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### Connection to the central gas supply

Attach unit's gas inlet ② direct to the wall outlet of the central gas supply (3.3...7 bar) using the appropriate adaptor hose.

**(i) Note:** Please use the low pressure gas supply setting (see below), as otherwise the warning 'CO<sub>2</sub> cylinder empty' is displayed at the start of insufflation (see p. 11).

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### Gas supply modes

The ThermoFlator® distinguishes between two different supply modes:

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**Installation and  
operating instructions  
Gas supply modes**

• **High pressure mode**

The high pressure mode corresponds to factory settings. When connected to an empty supply bottle or a low pressure supply, the left red LED of the supply pressure bargraph ② flashes. When the insufflation is started a warning sounds.

When the device is connected to the SCB system there is a permanent warning message on the PC or video monitor. Despite these warnings the device is fully functional when connected to a central low pressure supply.

• **Low pressure mode**

In this mode the right green LED of the supply pressure bargraph ② is always lit. There is no ‘cylinder empty’ warning.

**i Note:** If the unit is programmed for this mode but is connected to a high pressure supply it automatically switches over to the high pressure mode. However, this switchover only applies until power is switched off. Powering up the unit will always recall the supply pressure mode which was programmed.

The supply mode is set using service program P1 as described in section ‘Service programs’ on page 30. The programmed mode is saved even when the unit is switched off. The unit is delivered with high pressure mode as the default.

**Operating the unit**

Switch on the unit at its main power switch ①. The digital displays will illuminate shortly. Six short audible beep tones will sound to indicate that the unit is ready for use.

**i Note:** Once switched on, the unit awaits confirmation or a change of the last set values used. This is signaled by flashing of the corresponding LEDs in the bar displays. The LED in the START/STOP button also flashes.

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### Set value confirmation

To confirm all four set values (two set values each for two operating modes), one-time activation of any +/- buttons or the START/STOP button ⑯ is enough.

When the user confirms the setpoints a confirmation beep sounds. The beep will be heard whenever a confirmation via a front panel key or the SCB system took place.

Then the LED in the START/STOP button goes out and the LEDs in the bargraph displays are continuously illuminated. Subsequent activation of the START/STOP button then starts insufflation.

**Note:** The intermittent flow mode is selected as the initial setting.

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
Open the valve on the CO<sub>2</sub> cylinder about one-half turn counterclockwise.

**Caution:** The CO<sub>2</sub> cylinder must be placed vertically while operating the instrument, otherwise the functioning of the instrument cannot be guaranteed.  
The connected CO<sub>2</sub> cylinder must be secured to prevent it from falling over.

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### High pressure mode

Check the bar graph indicating the remaining gas supply ②.

- Left (red) LED flashing: CO<sub>2</sub> cylinder empty or central CO<sub>2</sub> gas supply connected
- All 6 green LED's light: adequate gas supply remains
- Only 3 green LED's light: remaining gas supply < 15 l
- Right (red) LED light: overpressure condition or wrong type of gas/gas cylinder

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### Low pressure mode

In this mode the right green LED of the supply pressure bargraph ② is always lit.

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Place a sterile CO<sub>2</sub> gas filter on the insufflation connection ⑬.

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**Caution:** To avoid contamination of the unit as a result of the reverse flow of CO<sub>2</sub> or body fluid, a sterile CO<sub>2</sub> filter must be inserted between the insufflation tube and unit inlet.

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
Fit insufflation tube to CO<sub>2</sub> gas filter. Connect the other end of the insufflation tube to the VERESS needle.

**i Note:** The VERESS needle can also be connected parallel to the insufflation tube using adaptor 20432031 and special connecting tube 20400044 (see accessories).

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
## Heating up the CO<sub>2</sub> gas

The ThermoFlator® enables the CO<sub>2</sub> gas to be pre-heated to body temperature. The risk of hypothermia during the insufflation of sizeable quantities of CO<sub>2</sub> can thus be reduced.

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### Connecting the Optitherm® heating element

- ❶ **Note:** The ThermoFlator® must be switched off when the Optitherm® heating element is connected.

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Connect the Optitherm® heating element ❷ to socket ❸ of the ThermoFlator® using the accompanying connecting cord.

- ❶ **Note:** If the connecting cord is disconnected, the ThermoFlator® can be operated without heating.

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## Installation and operating instructions

**Caution:** After sterilization in the autoclave, or after use of the Optitherm® heating element, a cool-down phase of 15 minutes must be observed before switching on again. Otherwise the varying cool-down behaviour of the temperature sensors may cause the alarm to be triggered during operation of the unit.

Pressing START/STOP button ⑯ automatically activates the heating element. The CO<sub>2</sub> gas flowing through is heated to a temperature of approx. 37 °C. The unit adapts to various flow rates automatically.

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Once the heating element is activated, indicator ⑫ will light:

- green if the heating element is working. The CO<sub>2</sub> temperature is preheated from room temperature to approx. 37 °C (99 °F).
- red in the case of overtemperature or a fault. An acoustic warning signal will be heard and the heating element is switched off; gas insufflation will be continued.

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**Note:** The overtemperature alarm resets itself if the gas temperature falls again to values below 41°C. Additionally, it can be shut off by pulling out the plug for the Optitherm® Heating Element.

For safety reasons, the alarm that is always triggered when a defective Optitherm® heating element is recognized (same alarm picture) cannot be reset. In this case, the unit must be switched off and switched back on again, either with a functional Optitherm® heating element or with no heating element at all.

In any case, the correct insufflation of (cold) gas is always guaranteed in spite of a recently sounded alarm.

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## Selecting the insufflation mode

The ThermoFlator® is equipped with two different insufflation modes:

- Insufflation initialization for the creation of pneumoperitoneum.
- Operation mode for insufflation during the operation.

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### Insufflation initialization

This operating mode is specially for positioning the VERESS needle and introducing pneumoperitoneum. In the ‘Insufflation initialization’ operating mode a continuous gas flow of 1 l/min and a pressure setpoint value of 15 mmHg are preset. The values can be changed with the +/– buttons ② and ⑩. However, the changes are not stored on leaving initialization mode.

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For insufflation initialization, press initialization mode switch ② until the indicator in the switch lights up. The preset setpoints for pressure and flow are indicated on displays ④, ⑤, ⑧ and ⑨.

**i Note:** The patient symbol ⑥ in the ‘Patient pressure’ display is extinguished in this operating mode.

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**i Note:** The automatic safe control identifies an unusually rapid increase in pressure (>2.5 mmHg/0.1 sec) as an obstruction in the insufflation line (e.g. needle tip in tissue instead of the abdominal cavity, kinked insufflation tube) and draws the user's attention to the error by an audible error signal as well as a flashing digital pressure display ⑤.

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- ❶ **Note:** As there is a slight underpressure in the peritoneum under physiological conditions, the unit may trigger an underpressure alarm when examination begins.
- ❶ **Note:** This operating mode also allows so-called low-pressure insufflation with an inflation pressure of up to 30 mmHg.

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### Operation mode

To enter operation mode, press initialization mode switch ② until the pilot lamp in the switch goes out.

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In operation mode, the patient pressure required can be preselected using +/- buttons ②.

The patient pressure preselected can be read on both the bar display ④ and the digital display ⑤. The flashing decimal point on the digital display indicates that the unit is in setpoint input mode.

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- ❶ **Note:** The setpoints in the 2 flow modes remain stored if the instrument is switched off or is switched over to a different mode.

Once it is switched back on, the unit awaits confirmation or change of the last setpoint values used (see Section ‘Set value confirmation’, page 18). Pushing a button (any +/- button or the START/STOP button) confirms the setpoint values.

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## Preselecting the gas flow

The ThermoFlator® is equipped with two different types of flow:

- Semicontinuous Flow
- Intermittent Flow

The unit always starts in the 'Intermittent flow mode' as standard.

Press on the ⑯ button bu to change between the two flow types.

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## Intermittent flow mode

In this operating mode, the gas flow is reduced to 0 during the pressure measuring phases.

It offers the advantage of precise measurement of intra-abdominal patient pressure.

To activate this operating mode, press the M button ⑯ until the 'Intermittent flow mode' indicator ⑰ and patient pressure symbol ⑥ light up.

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Then enter the required gas flow using ± buttons ⑮. The preselected value can be read on both the bar display ⑧ and the digital display ⑨. The flashing decimal point on the digital display indicates that the unit is in setpoint input mode.

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**(i) Note:** Naturally, the pneumoperitoneum can also be established in the "intermittent flow" operating mode.

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### Semicontinuous flow mode

In semicontinuous flow mode, the flow during the measuring phase is reduced not to 0, but only to 20 % of the setpoint set. This operating mode is particularly suitable for use with lasers and HF units since the semicontinuous flow largely ensures the absence of smoke.

**Note:** Complete discontinuation of gas flow triggers an alarm in this operating mode (audible error signal, flashing flow mode indication ⑨). The constant exchange of gas (at least 1 l/min) must therefore be ensured by the continuous removal of smoke by suction or by gas leakage. This facility prevents the line to the patient from being interrupted unnoticed and pressure monitoring from being distorted.

**Note:** 'Semicontinuous flow' is displayed when indicator ⑪ lights up. The patient symbol ⑥ in the 'Patient pressure' display is extinguished in this operating mode.

Then enter the gas flow required using the ± buttons ⑯.

The preselected flow can be read on both the bar display ⑧ and the digital display ⑨. The flashing decimal point on the digital display indicates that the unit is in setpoint input mode.

**Note:** The patient symbol ⑥ in the 'Patient pressure' display is extinguished in this operating mode.

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### SECUVENT® blow-off valve

In addition to electronic overpressure monitoring, the unit is fitted with an electromagnetic blow-off valve. As soon as an overpressure has continued for several seconds (approx. 5), an acoustic warning sounds. The red indicator ⑪ will light and the overpressure is reduced via the blow-off valve. The blow-off valve is enabled each time the unit is switched on, but can be disabled by the user.

**Caution:** If the blow-off valve is enabled; it is imperative that a sterile CO<sub>2</sub> filter be inserted between insufflation connection ⑬ and the insufflation tube to prevent the risk of contamination of the unit due to the reverse flow of gas or body fluid.  
The CO<sub>2</sub> filter must be replaced after each use.

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### Disabling the blow-off valve

Press the M button ⑯ for approx. 2 sec.; the blow-off valve indicator ⑪ will light permanently. The blow-off valve is now disabled. In the case of overpressure, the user will, however, still be warned by an audible signal.

To enable the blow-off valve, press the M button ⑯ again for approx. 2 seconds.

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
### Test for proper operation

**Warning:** Test this equipment prior to each surgical procedure to ensure that it functions correctly.

Preselect insufflation initialization mode.

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Quickly close the VERESS needle, activating monitoring of the pressure increase rate. An audible warning signal must be sounded and the pressure display ⑤ must flash.

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### Testing the blow-off valve

Leave the initialization mode by pressing button ⑯ and select ‘intermittent flow mode’.

Select the following pressure and flow settings: flow  $\geq 10 \text{ l/min}$ , pressure 15 mmHg.

Start insufflation and close the VERESS needle. If the patient pressure setpoint is exceeded for more than 5 mmHg for more than 3 seconds, a continuous warning signal must be sounded. After approx. 5 seconds the blow-off valve must release gas and the actual value of the patient pressure must now fall slowly until it is equal to the setpoint.

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### Initiating CO<sub>2</sub> insufflation

Set the volume display ⑩ to 000 by pressing the RESET button ⑯.

Introduce the VERESS needle and open the lever on its fitting. Press the START/STOP button ⑯ to initiate insufflation.

The insufflation initialization mode is recommended for creation of the pneumoperitoneum. The unit should then be switched over to operation mode.

**i Note:** HiCap® instruments which have an appropriate connection must be employed to use the maximum flow of the unit. The LUER-Lock connection fitted to the end of the insufflation tube must be removed in this case.

The insufflation pressure, flow and gas volume can be read continuously on displays ③, ⑥, ⑦, ⑨ and ⑩.

If the patient pressure rises above the preselected value, an audible warning signal (continuous tone) will be sounded. If the blow-off valve is enabled, the overpressure is actively reduced as well.

An audible warning signal (intermittent tone) is also sounded in the case of underpressure.

Close the valve on the CO<sub>2</sub> cylinder at the end of the procedure.

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### Service programs P1, P2

- Selecting the service programs

To select the service programs the ‘Volume Reset’ button ⑯ has to be pushed during power-up.

For about one second **P1** (Program 1) is displayed. Then, with **Hi SuP** (High Supply Pressure) or **Lo SuP** (Low Supply Pressure) the current mode of gas supply is displayed.

- Changing mode (service program P1)

This mode (Hi SuP <-> Lo SuP) can be changed with the + and – keys ⑯ for gas flow.

- Changing program (P1 <-> P2)

Pushing the + key for the pressure ⑯ selects the next service program **P2** (LED test program).

For about one second **P2** (Program 2) is displayed. This is the LED test program which runs in an endless loop.

The – key for pressure ⑯ selects the previous program **P1**.

- Saving mode or leaving service programs

To do this, switch off the unit. The user's inputs will be permanently stored.

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### Identical Devices on SCB

As from the R-UI Update CD 200900 01S243, the SCB R-UI system allows two or more identical devices of the same type to be operated in the SCB network and centrally controlled via the SCB software.

Up to 3 devices can be controlled via the SCB software. To be able to distinguish between the devices, they are each assigned an individual distinctive number.

Please attach the corresponding sticker (“#1”, “#2”, or “#3”) next to the SCB symbol on the front panel after you have assigned the relevant number.

You will find detailed information on this subject in section “Identical Devices on SCB” of the “KARL STORZ-SCB control NEO system” instruction manual.

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### Fuse replacement

Switch off the unit and remove the power plug from the electrical outlet.

Remove the power fuse holder ⑨ with a screwdriver or other tool.

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**Caution:** Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

	264320 20-1 100...240 V~
Power fuse	2 x T2AL250V

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Refit the power fuse holder ⑨.  
Reconnect the power cord.  
Test the unit for proper operation.

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**Warning:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" and the accompanying documentation.

**Warning:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

**Warning:** National laws and regulations must be observed.

**Warning:** Product damage: The procedures specified in the instruction manual have been validated by KARL STORZ. Only use these procedures.

**Warning:** Product damage: The reprocessing procedure must be validated by the user on site.

**Note:** The manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" can be downloaded or requested by visiting [www.karlstorz.com](http://www.karlstorz.com).

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**Optitherm® heating element**

For manual/thermal machine cleaning/disinfection (93 °C), use chemicals that have been specially approved by KARL STORZ for the reprocessing of medical devices. Use a soft brush for manual cleaning.

**Caution:** Do not damage heating wires!

**Caution:** When preparing and using the chemicals, follow the manufacturer's instructions paying close attention to proper dilution and working times. Prolonged immersion can lead to material changes.

The Optitherm® heating element can also be cleaned in an ultrasonic bath if there is heavier soiling. However the cleaning time should not exceed 3 minutes.

Sterilize the Optitherm® heating element at 134°C, preferably with the fractionated prevacuum procedure.

**Note:** The Optitherm® heating element absolutely must have cooled down to room temperature before use. If this is not the case, the error message 'Temperature exceeded' will appear and the heating cannot be started.

**Caution:** Do not remove the connector for cleaning, disinfection and sterilization.

**Caution:** Check for moisture residue in the contact connector before putting into operation. Dry with a lint-free disposable cloth if necessary.

**Inspection procedures for reusable silicone tubing sets before each sterilization**

**Warning:** For all cleaning and disinfection procedures on contaminated instruments, the local guidelines of the Employer's Liability Insurance Association or equivalent organizations must be observed.

**Warning:** A visual and functional check of the silicone tubing set must be carried out!

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- Visually inspect the silicone tubing set prior to each use to ensure that it is suitable for reuse. There should be no holes, cuts, dents or scratches in the silicone tubing. If the tubing set is damaged as described above, it must be replaced. If the tube set is discolored as a result of the reprocessing procedure, it must be removed immediately from the sterile supply cycle.

**Note:** Discolorations impede visual inspections for possible contamination.

- The integrity of the tubing set can be tested by attaching a syringe to one end of the tubing and doubling over the other end. Place the tubing in water. Force air into the tubing set using the syringe. If any bubbles appear in the water while forcing air through the tubing, the tubing set must not be reused.

**Note:** The leakage tester 13242 XL may also be used to check for leaks (see page 37).

**Cleaning, disinfection and sterilization of reusable tubing**

**Warning:** Risk of infection/product damage: The procedures specified in the instruction manual have been validated by KARL STORZ. Only use these procedures. The reprocessing procedure must be validated by the operator on site.

**Warning:** Reusable tubing must be cleaned, disinfected and sterilized before use and subsequent reuse using validated reprocessing procedures.

**Warning:** Prior to the sterilization of reusable tubing, particularly in cases of manual reprocessing, it must be ensured that all chemical residue has been rinsed off or removed (see also reprocessing steps).

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**Warning:** Disposable tubing sets, locks and hydrophobic bacteria filters must not be reused.

**Note:** For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the reprocessing of medical devices.

The list of KARL STORZ approved materials can be found on the Internet at ‘[www.karlstorz.com](http://www.karlstorz.com)’.

### Preliminary cleaning of reusable tubing

Place the reusable tubing in a container filled with a cleaning solution (as per manufacturer's instructions) immediately after use to prevent contaminants from drying on the surface of the tubing set.

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### Water Quality Requirements

Tap water can be used to prepare the cleaning and disinfectant solutions. When rinsing for the last time for example to remove any chemical residue, microbiologically pure/sterile water must be used.

### Reprocessing

#### Manual reprocessing of reusable tubing

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## Maintenance

### Cleaning, disinfection and sterilization

10. After reprocessing, inspect the tubes for residues and damage. If contaminants can still be detected, the reprocessing procedure must be repeated.

### Machine reprocessing of reusable tubing

1. Prior to machine reprocessing, carry out steps 1-5 of manual reprocessing (see page 36).

Thermal disinfection is preferred. The relevant national requirements and the  $A_0$  value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out, must take place in consultation with the manufacturer of the device.

**i NOTE:** If necessary, the instrument must be dried off afterwards by hand.

### Connection

2. In order to guarantee effective machine cleaning and disinfection, connect the tubing to the washer and disinfectant using the LUER connectors provided.
3. Attach the tubing to rinsing lances in order to ensure that it is rinsed out completely (with the plastic connectors, with the tube adapter end to the top of the rinsing lance).

### Checking the tubing set for leaks

**Warning:** Check the entire tubing set for leaks after reprocessing. In doing so, it is possible to detect a leak early, through which gas could escape.

To check for leaks use leakage tester 13242 XL, for example. To do so, connect one end to the leakage tester and block off the other end, for example, with your finger. If there is a leak in the tubing set, a pressure drop will be observed on the pressure gauge of the leakage tester.

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## Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual ‘Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments’. The method is to be selected taking into account the specific national requirements and in consultation with the device manufacturers.

### Steam sterilization in the fractionated prevacuum procedure

The medical device must be sterilized in its fully assembled state using a fractionated prevacuum procedure (EN ISO 17665-1) at 134 °C–137 °C with a minimum exposure time of 4 minutes and maximum of 18 minutes. These procedures are only suitable for thermostable instruments.

 **Note:** Use the fractionated prevacuum procedure (acc. to EN ISO 17665-1) observing all nationally applicable regulations.

 **Warning:** Successful sterilization is only possible on clean, disinfected surfaces.

 **Warning:** The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and validated.

 **Warning:** Any deviations from the recommended parameters for sterilization must be validated by the user.

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## Reprocessing limits

The end of the useful service life of the reusable tubing is largely determined by wear, reprocessing methods (chemicals) and damage resulting from use.

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## Maintenance and safety check

### Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

### Safety checks/repeat inspections as per IEC 62353

**Warning:** Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, with this device safety checks/repeat inspections as defined by IEC 62353 must be performed and recorded by a qualified electrician at least once a year.

### Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

### Electric measurements

- Inspect the device safety fuses
- Measure protective ground resistance in accordance with IEC 62353: The limit values can be found in the current standard.
- Measure earth leakage current in accordance with IEC 62353: The limit values can be found in the current standard.
- Touch current measured according to IEC 62353: The limit values can be found in the current standard.
- Measure patient leakage current in accordance with IEC 62353: The limit value for BF can be found in the current standard.

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### Test for proper functioning

Carry out the functional check as per the instruction manual (see section “Test for proper operation”).

### Documentation

The safety check/repeat inspection and results must be documented.

**Warning:** If any defects or shortcomings which could endanger patients, personnel or others are ascertained during these safety checks/repeat inspections, the device must not be operated until such time as these defects or shortcomings have been eliminated by qualified technical servicing.

**Note:** For detailed information on the range and implementation of the safety check/repeat inspection, please refer to the latest issue of the relevant Service Manual.

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### Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

### Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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### Repair program

For fiberscopes and equipment, individual repair is necessary. Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to:

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In other countries please contact your local KARL STORZ branch or authorized dealer.

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### Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff). We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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**Limitation of liability**

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations in the room where the unit is connected and operated comply with the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

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**Manufacturer's warranty**

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see "Subsidiaries" page), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

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### Troubleshooting

**Warning:** Always unplug the unit before all maintenance work!

#### Symptom:

Complete failure of the unit.

#### Possible causes:

No power from the power line.

Defective power fuse.

Temperature limiter is activated.

#### Remedy:

Check that there is electricity to the wall outlet.

Change fuses as described in the instruction manual. Use only fuses of the correct rating.

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No gas flow.

#### Possible causes:

CO<sub>2</sub> cylinder empty or CO<sub>2</sub> cylinder valve not fully opened. Insufflation deactivated.

#### Remedy:

Check the gas supply remaining on the indicator and replace the CO<sub>2</sub> cylinder if empty, or fully open the valve on the CO<sub>2</sub> cylinder. Press the "START/STOP" button to activate insufflation.

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#### Possible causes:

Leakage in the system of gas lines.

Failure in gas pressure regulation electronics.

#### Remedy:

Check system of gas lines, particularly all fittings, for leakage and replace, if necessary.

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#### Abhilfe:

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#### Symptom:

Heating element indicator lights red.

Acoustic warning signal sounds.

#### Possible causes:

Failure in the heating element.

Failure in heating regulation electronics.

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- Type of protection against electric shocks:  
Protection Class I
- Degree of protection against electric shocks:  
Applied part of type BF 

**According to IEC 60601-1-2:**

Please read the Electromagnetic Compatibility  
Information in the appendix (p. 49-61).

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## Technical description

### Directive compliance (for 264320 20-1)

#### According to Medical Device Directive (MDD):

This medical product belongs to Class IIb

This medical product bears the CE mark in accordance with MDD 93/42/EEC.

**i Note:** The code number after the CE mark indicates the responsible notified body.

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## Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated as repairable by their respective manufacturers.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

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
We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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**Funktionsschema Insufflation**

**Insufflation operational  
schematic**

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**Spare parts,  
recommended accessories**

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**Spare parts/accessories**

Item	Cat. no.
<b>Silicone tubing set</b> , sterilizable	<b>20400043</b>
<b>High pressure line</b> , american connector/ German standard connector, length 55 cm	<b>20400021</b>
<b>Same</b> , length 102 cm	<b>20400027</b>
<b>High pressure line</b> , american connector/ PIN-Index connector, length 55 cm	<b>20400022</b>
<b>Same</b> , length 102 cm	<b>20400028</b>
<b>High pressure line</b> , american connector/ american connector, length 102 cm	<b>20400025</b>
<b>Same</b> , length 55 cm	<b>20400029</b>
<b>Low pressure tube</b> , for the central CO <sub>2</sub> gas supply, length 300 cm	<b>20400026</b>
<b>High pressure line</b> , American connector/ ISO connector, length 55 cm	<b>20400221</b>
<b>Same</b> , length 102 cm	<b>20400222</b>
<b>CO<sub>2</sub> bottle</b> 1 l, empty, with German standard connector	<b>26400090</b>
<b>CO<sub>2</sub> bottle</b> 1 l, empty, with pin-index connector	<b>26400091</b>
<b>Open-end wrench</b>	<b>20400030</b>
<b>Optitherm® heating element</b>	<b>20432030</b>
<b>Power fuses T2AL250V</b> , pack of 10	<b>2027590</b>
<b>SCB connecting cord</b> , length 100 cm	<b>20090170</b>
<b>Power cord</b> (with grounding contact)	400 A
"Hospital Grade" <b>Power cord</b> (USA)	400 B
<b>Instruction manual</b>	96116011D
<b>Operating instruction quick reference guide</b>	96116011X

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**Accessories**

Item	Cat. no.
<b>Filter with insufflation tube,</b> for single use, sterile, pack of 10	031222-10
<b>Gas Filter,</b> for single use, sterile, pack of 25	031122-25
<b>VERESS Pneumoperitoneum Needle</b> with spring-loaded blunt stylet. LUER-Lock, Length: 10 cm	26120 J
<b>Same,</b> Length: 12 cm	26120 JL
<b>Same,</b> Length: 7 cm	26120 JK
<b>Same,</b> but 15 cm long, for setting up pneumoperitonea employing DOUGLAS-puncture in the case of adipose patients.	26120 JF
<b>Localization Needle,</b> LUER-Lock, O. D.: 0,8 mm, length 12 cm, for use with syringe 26003, 10 cc	26120 L
<b>Syringe,</b> 10 cc, LUER-Lock.	26003
<b>Spare barrel,</b> 10 cc., for 26003	26003 E
<b>Adaptor</b> for secondary connection of VERESS-needle	20432031
<b>HiCap® Trocar,</b> 11 mm	30103 LP
<b>High-pressure in-line filter</b>	20400032
<b>Leakage tester</b>	13242 XL

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## Appendix

### Electromagnetic Compatibility (EMC) Information

**Warning:** Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The Thermoflator® model 26 432020-1 corresponds to IEC 60601-1-2 :2007 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The Thermoflator® model 26 432020-1 is a Group 1 unit (as per CISPR 11). Group 1 contains all the ‘equipment and systems which generate or use RF energy only for their internal functioning’.

**Note:** The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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## Appendix

### Electromagnetic Compatibility (EMC) Information

**Warning:** The ThermoFlator® model 26432020-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

**Warning:** The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

**Warning:** The use of accessories, transducers and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the ThermoFlator® model 26432020-1. The accessories, transducers and cables listed below have been shown to comply with the requirements of IEC 60601-1-2. When using accessories, transducers and cables other than those specified here, it is the responsibility of the user to ensure that they comply with IEC 60601-1-2.

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**Accessories, transducers and cables which have been shown to comply with IEC 60601-1-2:**

Type	Shielded	Length [m]	Ferrite	Used for
PA	No	> 3	No	Potential equalization
204320 30	Yes	2.5	No	Optitherm® model 204320 30
200900 30	Yes	4	No	SCBcom model 200900 30
Power cord	No	3	No	Power supply

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**Table 1**  
**Guidance and manufacturer's declaration – electromagnetic emissions**

Thermoflator® model **26 4320 20-1** is intended for use in the electromagnetic environment specified below.  
The customer or user of the Thermoflator® model **26 4320 20-1** should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Thermoflator® model <b>26 4320 20-1</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Thermoflator® model <b>26 4320 20-1</b> is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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**Table 2**  
**Guidance and manufacturer's declaration – electromagnetic immunity**

Thermoflator® model 264320 20-1 is intended for use in the electromagnetic environment specified below. The customer or user of the Thermoflator® model 264320 20-1 should ensure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % $U_T^*$ (>95 % dip in $U_T$ ) for 0.5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	Complies <5 % $U_T^*$ (>95 % dip in $U_T$ ) for 0.5 cycle  Complies 40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  Complies 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  Complies <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment.  If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

\* Note:  $U_T$  is the a.c. mains voltage prior to application of the test level.

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**Table 4**  
**Guidance and manufacturer's declaration – electromagnetic immunity –  
for equipment and systems that are not life-supporting**

Thermoflator® model 26 4320 20-1 is intended for use in the electromagnetic environment specified below.  
The customer or user of the Thermoflator® model 26 4320 20-1 should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3 V <sub>rms</sub>	Portable and mobile RF communications equipment should be used no closer to any part of the Thermoflator® model 26 4320 20-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5 / 3] \sqrt{P}$ 150 kHz to 80 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = [3.5 / 3] \sqrt{P}$ 80 MHz to 800 MHz $d = [7 / 3] \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> . Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects and persons.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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**Table 6**

**Recommended separation distances between portable and mobile RF communications equipment and the  
Thermoflator® model 264320 20-1**

The Thermoflator® model 264320 20-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters whose maximum nominal power is not listed in the above table, the recommended separation distance d in meters [m] can be determined using the equation belonging to the relevant column, where P is the maximum nominal power of the transmitter in Watts [W] as specified by the manufacturer of the transmitter.

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KARL STORZ Endoscopy Canada Ltd.  
7171 Millcreek Drive, Mississauga  
Ontario, L5N 3R3, Canada  
Phone: +1 905 816-4500, Fax: +1 905 858-4599  
Toll free phone: 1-800-268-4880 (Canada only)  
Toll free fax: 1-800-482-4198 (Canada only)  
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100, Fax: +1 424 218-8526  
Toll free phone: 800 421-0837 (USA only)  
Toll free fax: 800 321-1304 (USA only)  
E-Mail: info@karlstorz.com

KARL STORZ Veterinary Endoscopy America, Inc.  
175 Cremona Drive  
Goleta, CA 93117, USA  
Phone: +1 805 968-7776, Fax: +1 805 685-2588  
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.  
815 N. W. 57th Avenue, Suite 480  
Miami, FL 33126-2042, USA  
Phone: +1 305 262-8980, Fax: +1 305 262-8986  
E-Mail: info@ksela.com

KARL STORZ Endoscopia  
Miramar Trade Center,  
Edificio Jerusalén, Oficina 108,  
La Habana, Cuba  
Phone: +53 72041097, Fax: +53 72041098

KARL STORZ Endoscopia México S.A. de C.V  
Lago Constanza No 226 , Col. Chapultepec Morales  
D.F.C.P. 11520, México, México  
Phone: +52 55 1101-1520, Fax: +52 55 55450174

KARL STORZ Marketing América Do Sul Ltda.  
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,  
CEP 04534-011 São Paulo, Brasil  
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680  
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.  
Zufriategui 627 6º Piso, B1638 CAA – Vicente Lopez  
Provincia de Buenos Aires, Argentina  
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773  
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopie Norge AS  
Østensjøveien 15B  
0661 Oslo, Norway  
Phone: +47 6380 5600, Fax: +47 6380 5601  
post@karlstorz.no

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KARL STORZ ENDOSKOPE  
East Mediterranean and Gulf S.A.L.  
Block M, 3rd Floor  
Beirut Souks, Weygand Street  
2012 3301 Beirut, Lebanon  
Phone: +961 1 999390, Fax: +961 1 999391  
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.  
P.O. 6061, Roggebaai 8012, Cape Town, South Africa  
Phone: +27 21 417 2600, Fax: +27 21 421 5103  
E-Mail: info@karlstorz.co.za

## Sociedades distribuidoras

KARL STORZ Endoskopie (Shanghai) Ltd.  
Chengdu Branch  
Room 805, 8F Jin Jiang International Building  
No. 1 West Linjiang Road, Wuhou District,  
Chengdu, 6100414, People's Republic of China  
Phone: +86 28 86587977, Fax: +86 28 86587975  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shenyang Branch,  
Room 2001-2005, 20F N-MEDIA International Center,  
No. 167 Youth Avenue, Shenhua District,  
Shenyang, 110014, People's Republic of China  
Phone: +86 24 23181118, Fax: +86 24 23181119  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Guangzhou Branch  
Room 02B, 35F Teem Tower,  
No. 208 Tianhe Road, Tianhe District,  
Guangzhou, 510620, People's Republic of China  
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.  
I-Park Bundang1, 102dong, 8F  
239 Jeongjail-Ro, Bundang-Gu  
Seongnam, 463-859, Korea  
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299  
E-Mail: stephen.Lee@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.  
12F, No. 192, Sec. 2, Chung Hsin Rd.,  
Sindian District, New Taipei City, Taiwan  
Phone: +886 933 014 160  
Fax: +886 2 8672 6399  
E-Mail: mingwang@karlstorz.tw

KARL STORZ Endoscopy Japan K. K.  
Stage Bldg. 8F, 2-7-2 Fujimi  
Chiyoda-ku, Tokyo 102-0071, Japan  
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633  
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

**KARL STORZ GmbH & Co. KG**

 Mittelstraße 8  
78532 Tuttlingen  
Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 (0)7461 708-0  
Telefax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
Web: [www.karlstorz.com](http://www.karlstorz.com)