



COMPARE TO		
<i>drive</i>	BREWER®	CLINTON®
13035	43460	M21

## Mayo-Instrument Stand

### Single Post

Item #	UOM
13035	1/cs

- Single post stand with “U” base design provides support and versatility.
- Removable stainless steel tray measures 19" x 16.5".
- Tray height adjusts from 32"– 49" with lock.
- Two, 3" casters provide additional mobility.

specifications	
Tray Dimensions	19" (W) x 16.5" (D)
Height	32" - 49"
Base Dimensions	22" (W) x 18.5" (D)
Weight	9 lbs.
Carton Shipping Weight	11 lbs.
Warranty	Limited Lifetime
accessories	
13035T	Replacement Tray 1/cs